

A Miniature Multisensor Biosignal Data Recorder and its Evaluation for Unsupervised Parkinson's Disease Data Collection

Chris Bailey, Garry Hollier,
 Anthony Moulds, Michael Freeman
 Department of Computer Science
 University of York, UK
 chrisb@cs.york.ac.uk
 anthony.moulds@york.ac.uk
 garry.hollier@york.ac.uk
 mike.freeman@york.ac.uk

Jim Austin,
 Alex Fargus, Thomas Lampert
 Cybula Ltd,
 York, UK
 austin@cybula.com
 alex.fargus@york.ac.uk
 tom.lampert@hotmail.com

Abstract—An improved miniature biosignal data sensor and recorder device is described, (NAT-1-4G) with 3-axis accelerometer, and a 2K Sa/sec all-channel recording capacity of 24 hours or more with a single zinc-air battery cell. Like the previous NAT-1 prototype device, this measures less than 18x22x10mm and weighs less than 2.3 grams, including the battery. In this paper we describe the device in detail, and introduce the presentation of tremor data measurement captured in the context of Parkinson's disease fore-arm monitoring. The NAT-1-4G device has already achieved translation to commercialization and we expect it to become commercially available in the near future.

Keywords- *Neurophysiology; Bio-signal sensors; Medical sensors, Parkinson's Disease.*

I. INTRODUCTION

It is well understood that use of biosignal data acquisition is increasingly important in many application scenarios, not the least of which are biomedical applications. Often, such measurements are taken in wired or wireless umbilical modes, in other words, within a clinical evaluation setting, with data captured and analyzed over relatively short time windows, and in unnatural settings.

However, the ability to perform ambulatory monitoring of patients provides the possibility of long-duration data capture of bio-parameters in a normal living situation or work-place. This has been an aim for many decades, and has developed from early magnetic tape based data capture [1], to digital systems [2], custom integrated circuits [3][4] and more advanced medical data recorders [5][6]. Such capabilities are identified by many clinical researchers as being desirable. The opportunity to learn more about medical conditions as well as the condition of individual patients themselves is seen as a major motivator for developing suitable devices. This is very true of application in the domain of Parkinson's Disease and similar tremor-related medical conditions, where often the primary mode of data capture is limited to a supervised scenario [7][8][9]. The importance of gathering continuous data for drug management, establishing and gauging long-term prognosis, and the personal reassurance for patients afforded by suitably

presented feedback from monitoring systems, cannot be underestimated.

The "NAT" (Neural Activity Tracker) project aims to produce a multi-purpose data sensing and recording solution that is extremely small, lightweight and having a recording capacity of days to weeks, dependent upon the selection of parameters such as sample frequency. The NAT-1 device [10] introduced our first solution to this problem in 2013. Since then we have enhanced and improved our design, resulting in NAT version 1-4G. This newer device includes multiple improvements, including the ability to alter the signal sensitivity of individual recording channels to adapt to a given bio-signal type, including ECG/EEG/EMG/EOG (Electro-encephalography (Electro-myography, Electro-cardiography, Electro-oculography), and also auxiliary sensor modules (such as temperature for instance). The NAT-1-4G also retains its capability for tri-axial accelerometer data capture, and in an accelerometer-only recording mode allows 500 samples per second across three simultaneous axes, with recording times up to 3 days in current prototypes.

This paper presents the NAT-1-4G in technical detail in Sections II and III, whilst Section IV documents some initial experimental applications of the device with a Parkinson's Disease tremor collection scenario. Note that the intention of this paper is not clinical, we do not make any clinically definitive claims about the data collected, but will observe and determine the suitability of the device for capturing such data, and for showing characteristics typically symptomatic of a Parkinson's Disease subject.


Section V highlights the additional resources developed alongside NAT, including docking station, wrist-mount, and the infra-red data signature time-stamp daughterboard. A brief state-of-the-art is given in Section VI, and conclusions are presented in Section VII.

II. A SMALL FORM-FACTOR DATA-RECORDER

The NAT is a device family of only 18x22mm in size, less than 10mm in height (when using a zinc-air disposable power cell). It weighs less than 2.3 grams, which means it is attractive in applications where regulatory constraints apply

(e.g., use with small animals). It is so light in weight that it quite unobtrusive as a human-subject wearable body-sensor. This low-weight attribute also means that multiple devices can be worn individually or in small groups where appropriate, without causing encumbrance of the subject's normal movements and behavior. We have verified this in a fashion, with team members wearing the device continuously for up to three days (using a 3D-printed wrist-capsule) with no issues.

A NAT device is shown in the image of Figure 1, with a ball-point pen of normal size for scale. In the photograph provided in Figure 1, we can see some interesting features. The major part of the device has a profile of less than 4 mm. Also visible in Figure 1 is a specially designed battery clip for housing a zinc-air cell (as used in hearing aids). This is the gold-plated metal structure. The use of flatter button cells can be envisaged to make a smaller profile possible, though this has weight implications. We also note the potential for ultra-thin flexible lithium polymer rechargeable batteries, which have a similar footprint and only 2mm height profile,



that would align well with the NAT PCB module in situ.

Figure 1. Top View of NAT Device




Figure 2. NAT Device – Sub-component Block Diagram.

Near the tip of the ball-point pen, one can see a small low-profile connector socket, which is actually a mezzanine daughterboard connector port. This permits a range of possible extension modules. Figure 2 shows the system level block diagram. As for our previous NAT-1 device [10], this device comprises of three key components – a flash memory chip (largest chip in Figure 1), a proprietary CPU (mid-sized chip of Figure 1), and a proprietary MEMS (Micro-Electro-Mechanical System) device for accelerometer (smallest chip of Figure 1). Additional analogue front-end components provide appropriate signal conditioning for the signal ranges typically encountered in a range of biosensor and biomedical applications. This is augmented by the introduction of programmable gain on a per-channel basis, with signal range up to +/-8mV, and ability to scale to signals in the sub-1mV range without appreciable noise.

III. NAT-1-4G DEVICE SPECIFICATIONS

The NAT-1-4G has the specifications as outlined in Table I. The device has a wide range of possible sample rates, ranging from 100 Sa/sec to 2 KSa/sec via the user interface software application. At 2 KSa/sec, the device consumes 4.8mA of current from a single 1.4-volt cell, and can record for up to 12 hours. An 8-Gbit flash option is possible (NAT-1-8G) and would have up to 24 Hours of recording capacity at maximum rates. At lower sample rates, the capacity of the flash is extended to many days, for example 6 days at 100 Sa/sec for accelerometer only. This compares well with reported state of the art examples [3,4], given that the system is comprised of readily available commodity integrated circuits. The device has three important connection mechanisms, these being the analog input connector (angled connector block to left of device in Figure 3a), the daughterboard extension socket (Figure3b) and the docking-shoe connectors (seen in Figure 3a)

A particular feature of NAT-1 was the use of a single zinc-air power cell. Such cells are widely used in hearing aids, and have an active power-delivery life-span of around two weeks, after which the cell begins to lose effectiveness. We have retained this mode of power provision in NAT-1-4G. However the possibility of using alternative power cells is being actively investigated at present. To date, the Zinc-Air modules appear to give the best power density for size and weight with considerable record times up to 2 weeks.


IV. EVALUATION METHODOLOGY

Previously, most of the testing of the prototype device has been limited to test scenarios including (a) EEG and Accelerometer data collection from mice in live test scenarios undertaken by researchers at the University of Aberdeen, (b) a preliminary EMG capture evaluation at the University of York, and a (c) high-G test setting using a golf-club attachment to capture golf swing behaviors in terms of club rotation, side-movement and swing-path motion [10]. In all cases, the device was used untethered and powered exclusively with a single zinc-air cell. This


work was aimed at gaining an initial scope of capabilities of the device. In this paper, with the newer NAT-1-4G, we describe the use of the device in capturing tremor motion in a Parkinson’s disease subject, and we compare this to data captured in the same way for healthy control subjects. This is a limited test-scenario owing more to engineering evaluation than clinical study, and should not be considered as a practical clinical evaluation study.

The experimental method consisted of a wrist-mounted device (see later section VI for details and photographs) which in most cases is worn continuously for several days, and only being removed where not appropriate for the device operation (e.g. whilst showering). All three accelerometer channels were sampled at 500 Sa/Sec with a resolution of 8 bits per channel for a full-scale range of +/- 2000 milli-G. An important aspect of the data collection is that there is no knowledge of activities during data collection (it is unsupervised).

Unlike short laboratory assessed (supervised) tremor measurements, where a series of directed or constrained actions are monitored for a period of perhaps 1 hour, here the data was collected blindly for up to 72 hours. This presents particular issues for data interpretation which we discuss further in later sections. The purpose of the evaluation was to make a straightforward determination that the NAT-1-4G device is (or is not) capable of gathering potentially clinically-useful tremor data during long unsupervised periods of use.



(a)



(b)

Figure 3. (a) Overhead View of NAT and Daughterboard Connector, and (b) Analogue Connector Port on NAT Device Underside

In order to understand our evaluation, it is useful to first of all exemplify the typical data characteristics likely to be observed in such an experimental study. A common aspect of Parkinson’s tremor is the presentation of an excess of frequency components in the 5 to 7 Hz range as compared to typical healthy subjects.

The example of Figure 4, reproduced from a recent publication [11], illustrates this quite well. It is seen that a control subject has a relatively broad range of frequency components in movement measurements, whilst a Parkinson’s patient exhibits a classical dominant frequency peak centered around the 5-6 Hz frequency, (often described as the unilateral or bilateral resting tremor) which is attributed to motor-neuron induced involuntary movements in the instrumented limb. These components usually have noticeably larger magnitudes than typical average motions generated by intended limb motion.

TABLE I. NAT-1-4G SPECIFICATIONS

Parameter	Limits	Units
Analogue inputs	4	channels
Bits per analogue channel recorded	11	bits
Mems Accelerometer	3	Axis
Bits per Accel. Axis	8	bits
Sample rate (max)	2000	x 4 ch
Max Current	2K Sa/s 500 Sa/s	4.8 2.4 mA mA
Data Capacity	4 or 8	Gbits
All-Channel Recording Time at 2K Sa/Sec	12 or 24 ^a	Hours
Analogue Range	± 500 min ± 8000 max	uV
(previous NAT)	± 1000 ^b	
Accelerometer range	Selectable 2 / 8	G (G-force)
Accelerometer sensitivity	18 / 72	mG (G-force)

^a Recording time for all channels. Revised device has accelerometer-only option with 6 day record time at 100 Sa/Sec.

^b Original NAT had fixed +/-1mv voltage range. Revised NAT range is channel programmable up to +/- 8mv.




Figure 4, Example of Parkinson’s Tremor components centered around 5 to 6 Hz for control (top) and affected patient (bottom).
 Reproduced from Meigal et al [11]

It is important perhaps to note that there is a significant difference between a resting limb exhibiting tremor (as might be observed under a controlled measurement scenario) and a limb measured over a significant timescale with unsupervised ‘every-day’ behaviors. The challenge for long-term unsupervised monitoring is to distill reliable indicators from inherently component-rich data sets. NAT-1-4G can of course be used for both a controlled resting limb study and an unsupervised data recording scenario.

In terms of existing devices, of which there are too many to present a comprehensive survey here, we have surveyed a subset of current devices that have succeeded in achieving commercial translation into common use, as tabulated in Table II, for comparative purposes. Several of these are wirelessly tethered devices, and several are non-wireless devices, including NAT-1-4G itself.

The Neurologger-2A device has some similar sets of capabilities to the NAT-1-4G, as might be expected as it is a similar onboard storage recorder, whereas the wireless devices have significant demands for power, requiring larger batteries and more circuitry. It can be seen that NAT-1-4G provides the smallest, lightest, and longest recording lengths for a device integrating both accelerometer and analogue sensor (e.g. EEG) channel recording as standard. Neurologger-2 with add-on accelerometer board has comparable weight and size to NAT-1-4G. The NAT device has a large signal range and is programmable to accept analogue signal inputs up to +/-8mV and can be reduced to < +/- 500uV for high-sensitivity .

Ambulatory monitoring has importance for the role of extended out-patient monitoring of EEG as a diagnostic aid and guide to treatment [12][13]. Recent work in this field includes many wireless approaches to ambulatory monitoring, but these are effectively 'tethered' via a secondary monitoring or recording station [14][15]. Recent literature reports an 8-channel EEG sensor node which “measures 35mm x 30mm x 5mm excluding Li-ion battery” [14], and power consumption in the range of 27-42mW., equivalent to around 3 to 5 mW per channel.

TABLE II. COMPARISON WITH SOME RECENTLY REPORTED DEVICES

Device & Type	SIZE & WEIGHT	Battery Type & Record Time	Sensors and Inputs	Sample Rates	Data Storage	Refs if any
MROBERTS HYBRID Bluetooth Tethered	87x45x17 mm 74 grams	Lithium Polymer 60 hrs at 100 Hz	2G tri-axial accelerometer triaxial Gyroscope	100 Sa/sec	SD card	[16]
Neurologger 2A Untethered EEG logger	22x15x5 mm 1.7 grams with battery and 3-axis accelerometer	2x ZN10 battery 124 hrs at 100 Hz 32 hrs at 400 Hz	4 channel analog input extra board provides triaxial accelerometer with additional weight (about 0.4 gram).	100 Sa/sec up to 19 kSa/sec for one channel 1000 Hz for accelerometer	On board flash Up to 1Gbit	[17]
G-Link Wireless tethered	58x43x21 mm 40 grams	Lithium polymer Record time not stated	2 G/10 G Triaxial accelerometer, temperature sensor	2 channel (accel/temp) 4 kSa/sec 4 channel 2 kSa/Sec	Onboard 2Mbyte	[18]
NAT-1-4G Untethered EEG/EMG/ECG and Accel Logger	22 x18x10 mm < 2.3 grams with battery	Zinc-Air hearing aid battery 9 hours at 2000 Sa/Sec 72 hours at 250 Sa/sec >1 week at 100 Hz	2 G/8 G triaxial Accelerometer 4 channel +/-8 mV analog in	Up to 2000 Sa/sec	Built in Flash 4 or 8 Gbit.	This paper

V. EVALUATION RESULTS

After gathering data from a total of four subjects, various data processing approaches were explored to derive useful data. The initial raw data sets contain a number of artifacts that need to be managed. Inevitably the accelerometer range encounters clipping (saturation), since the sensitivity of the accelerometer scale needs to be fine enough to detect tremor motions, even though other motions in unsupervised data collection can be more dynamic. A second issue is that of constant-G offsets caused by orientations of the limb (in this case the wrist) in different postures for instance. Both of these aspects may be observed in Figure 5a, which shows the raw data collected over a 48-hour period, and clearly exhibits variable constant-G offsets as well as saturation.

Figure 5a highlights a further aspect of interest, colored in red on the trace. In the data survey we have conducted we have observed that it is possible to automate identification of ‘quiet periods’ where activity is relatively low. These are analogous to a resting limb, though in fact these so-called quiet periods contain non-resting activity, but less so than other segments of the data.

By applying very-low-frequency component extraction and adjustments to the raw data, it is possible to remove the constant-G offsets from the data, whereupon it is possible to plot the magnitude of acceleration of the data on each channel without bias due to such offsets. This is illustrated by Figure 5b, which shows the normalized magnitude in red against the whole data set in pink. Figure 5c shows a small segment of the data with the raw data magnitude plot versus the offset adjusted case.

Once these treatments are applied to the gathered data, it becomes possible to generate power-spectra plots of magnitude versus frequency components using FFT analysis. We may then plot each accelerometer axis as well as a combined power spectra plot, which combines motion from all three axes. Figure 6a shows such a collection of plots for a healthy subject, whilst Figure 6b shows the same for a known Parkinson’s disease volunteer. It is already noticeable that the Parkinson’s disease plots show the classical excess of activity in the 5 to 6 Hz region.

Figures 7a and 7b show the combined data sets for four human volunteers for comparison. The subjects included three healthy individuals and one Parkinson’s disease diagnosed subject. Data was collected for each wrist (left/right) on separate occasions. When only the quiet periods of the data sets are analyzed, the components of the data become less ‘cluttered’. That is not to say there is not much to be discovered by analyzing a more complex data context – but rather that there is much to be done to develop useful analytical models from such data sets.

Figures 8a and 8b show examples of non-selective data analysis (i.e. all data from a multi-day data capture run) and using the so-called ‘quiet periods’ only.

In the healthy subjects the plots have a clearly defined roll-off toward 20 Hz and if the slope of the characteristic was subtracted from these plots there would be a fairly uniform response. In the case of the Parkinson’s subject however, one can observe an obvious excess of frequency

activity in the 5 to 7 Hz range and this correlates well with our expectation if the device were capable of effective collection of tremor data. We note that the left and right hand characteristics of the Parkinson’s disease volunteer are different. This may reflect the fact that the dominant limb and non-dominant limb have different use in every-day activity and/or a function of the tremor being non-symmetrically presented.

All of the subjects exhibited a strong signal in the area of 1 to 2 Hz, which we attribute to normal behavior – most likely arm swing due to walking and gait behavior. This will be subject to further investigation in future work planned by the group. As noted earlier this is not intended to be a clinical study – rather an observation of the suitability of our device for long term signal capture in this domain – can we capture data that is useful in differentiating tremor related micro-movements? The data captured seems to confirm the value of the device for such work.

VI. DOCKING STATION, ANALYSIS TOOLS, AND EXTENSION DAUGHTER BOARDS

The NAT device has the advantage of being provided with a bespoke docking station, for download of data from flash via USB. This is illustrated in Figure 9. Once downloaded, the data can be processed using Cybula Ltd Signal Data Explorer software suite, which can be trained to perform auto detection and classification of signal behaviors and events [10]. A final notable aspect of the current NAT-1-4G system suite is the infra-red time-code recording daughterboard. This extremely small add-on allows recording of an infrared pulse code stream along-side the analogue channels and accelerometer (see Figure 10 for photograph of IR device fitted on NAT main unit).

Figure 11 shows the wrist-mount prototype capsule in various views when worn. The wrist module is quite compact but could be further reduced in size. With certain power cells the format and size of the capsule could certainly be reduced considerably. These early prototypes were developed using 3D rapid prototyping (close up view shown in Figure 12). We are currently working on more sophisticated one-piece wrist-mount modules with integral strap, and using more flexible plastics.

Future developments for the daughterboard connector socket could include RF telemetry functions. It is also possible to use the daughterboard connector to act as a micro-backplane to permit multiple NAT devices to be ganged and synchronized to perform something like 16 or 32 channel acquisition, thus linking multiple NATs in a synchronized data capture mode. This is particularly interesting for multi-limb monitoring where relationships between respective limb data-sets might be of interest.




Figure 5. (a) Raw Data Collection over 48 hours. (where red sections represent 'quiet' periods). (b), Raw data and Magnitude after Constant-G Correction (where red represents adjusted data).




Figure 5c, Magnitude of Acceleration with Correction (blue - original data, green - offsets removed)




Figure 6a, Three-Axis Power Spectra (Non-Parkinson) (right - dominant - hand), showing channels 1,2, and 3, plus combined magnitude, linear vertical scale.






Figure 6b, Three-Axis Power Spectra (Parkinson) (right - dominant - hand, showing channels 1,2, and 3, plus combined magnitude, linear vertical scale)



(a)



(b)

Figure 7b, Power Spectra for Parkinson's Case and Controls: (a) Left Wrist, (b) Right Wrist. (Parkinson's subject red, blue and green are controls, right wrist Data Collection, log vertical scale)





Figure 8, Frequency vs. Magnitude for All , (a) without selectivity, and (b) with 'quiet period' selectivity (right wrist, trace colours as noted in Figure 7.)



Figure 9. NAT USB docking station




Figure 11 NAT Wrist-Mounted Early Prototype Capsule




Figure 10 NAT1-IRDB Infra-red Receiver Daughter-Board, Shown mounted on NAT device



Figure 12, Close-up of NAT Capsule Prototype.

VII. CONCLUSIONS AND ONGOING WORK

The improvements made to the NAT-1-4G device are undoubtedly beneficial for its wider use in application areas such as EEG, EMG and in more complex mixed monitoring scenarios. This includes combined tremor-motion and EMG capture in patients with conditions such as Parkinson's disease: a key area of interest at the present time. Research continues within our group on more sophisticated data analysis and algorithmic treatment of gathered data in long-duration unsupervised contexts. The use of low sample rates for accelerometer data recording of the order of 7 days is possible with the newer device, and we hope in future to evaluate the tradeoff between data quality and sample rates to establish a clinically useful long-duration recording configuration. The York team [19] and partners CYBULA Ltd. [20] are very interested in developing collaborative partnerships and would welcome enquiries from prospective evaluators and end-users in the clinical domain at the present time.

This work, forms the first phase of a longer-term research project in which we expect to continue to collect data in increasingly sophisticated scenarios with continuous data capture over months and possible years. The particular aim is to achieve clinically useful continuous unsupervised and untethered monitoring in every-day environments. The availability of rich data sets in this area from a single device will be extremely valuable and will also guide development of the next generation NAT devices.

ACKNOWLEDGMENTS

The NAT-1 and NAT-1-4G devices have been developed through a number of collaborative activities, including Cybula Ltd, University of York, and collaboration with the University of Aberdeen. We are also grateful for the support of the UK Government's Technology Strategy Board funding agency in supporting the development of these devices to date, and for funding from the University of York Strategic Investment Fund.

REFERENCES

- [1] J. R. Ives, and J. F. Woods. "4-channel 24 hour cassette recorder for long-term EEG monitoring of ambulatory patients." *Electroencephalography and clinical Neurophysiology* 39, no. 1 (1975): pp. 88-92.
- [2] A. L. Vyssotski et al. "Miniature neurologgers for flying pigeons: multichannel EEG and action and field potentials in combination with GPS recording." *Journal of neurophysiology* 95, no. 2 (2006): 1263-1273.
- [3] R. F. Yazicioglu, P. Merken, R. Puers, and C. Van Hoof, "A 200 μ W eight-channel acquisition ASIC for ambulatory EEG systems". *Solid-State Circuits Conference, 2008. ISSCC 2008, San Francisco. IEEE International Digest of Technical Papers, Feb 2008*, pp. 164,603.
- [4] N. Verma, A. Shoeb, J. V. Guttag, and A. P. Chandrakasan. "A micro-power EEG acquisition SoC with integrated seizure detection processor for continuous patient monitoring." In *2009 Symposium on VLSI Circuits, 2009*, pp. 62-63.
- [5] E. Waterhouse. "New horizons in ambulatory electroencephalography." *Engineering in Medicine and Biology Magazine, IEEE* 22, no. 3 (2003): pp. 74-80.
- [6] F. S. Jaw, Y. L. Tseng, and J. K. Jang. "Modular design of a long-term portable recorder for physiological signals." *Measurement* 43, no. 10 (2010): pp. 1363-1368.
- [7] W. Maetzler, J. Domingos, K. Srulijes, J.J. Ferreira, and B.R. Bloem. "Quantitative wearable sensors for objective assessment of Parkinson's disease." *Movement Disorders* 28, no. 12 (2013): pp. 1628-1637.
- [8] M. A. Hobert, W. Maetzler, K. Aminian and L. Chiari. "Technical and clinical view on ambulatory assessment in Parkinson's disease." *Acta Neurologica Scandinavica* (2014), Volume 130, Issue 3, pp. 139-147.
- [9] R. LeMoyne, T. Mastroianni, and W. Grundfest. "Wireless accelerometer configuration for monitoring Parkinson's disease hand tremor." *Advances in Parkinson's Disease 2* (2013): pp. 62-67.
- [10] C. Bailey et al. "A Miniaturized 4-Channel, 2Ks/sec Biosignal Data Recorder With 3-Axis Accelerometer and Infra-red Timestamp Function." In *SENSORCOMM 2013, The Seventh International Conference on Sensor Technologies and Applications, 2013*, pp. 213-219.
- [11] A. Meigal et al. "Linear and nonlinear tremor acceleration characteristics in patients with Parkinson's disease." *Physiological measurement* 33, no. 3 (2012): pp. 395-412.
- [12] S. A. Rizvi, J. F. Téllez-Zenteno, S. L. Crawford, and A. Wu. "Outpatient ambulatory EEG as an option for epilepsy surgery evaluation instead of inpatient EEG telemetry." *Epilepsy & Behavior Case Reports* 1 (2013): pp. 39-41.
- [13] H. J. Faulkner, H. Arima, and A. Mohamed. "The utility of prolonged outpatient ambulatory EEG." *Seizure* 21, no. 7 (2012): pp. 491-495.
- [14] S. Patki et al. "Wireless EEG system with real time impedance monitoring and active electrodes." In *IEEE Biomedical Circuits and Systems Conference (BioCAS), 2012*, pp. 108-111.
- [15] S. Mitra, J. Xu, A. Matsumoto, K. A. A. Makinwa, Chris Van Hoof, and R. F. Yazicioglu. "A 700 μ W 8-channel EEG/contact-impedance acquisition system for dry-electrodes." In *2012 Symposium on VLSI Circuits (VLSIC), IEEE, 2012*. pp. 68-69.
- [16] McRoberts Hybrid, commercial product data-sheet, <https://www.mcroberts.nl/files/Flyer%20Hybrid%20En2.pdf> f, url checked 21st July 2014.
- [17] Neurologger 2A Product Data Sheet, <http://www.evolocus.com/Products/Neurologger2A.pdf> url checked 21st July 2014.
- [18] G-Link-LXRS datasheet, Lord Corporation MicroStrain® SensingSystems (USA), www.microstrain.com url last checked 21st July 2014.
- [19] Advanced Computer Architectures Group (ACAG). www.cs.york.ac.uk/arch , last accessed 01/10/2014.
- [20] Cybula Limited. UK Limited Company. www.cybula.com , last accessed 01/10/2014.